



S-FRAME

1



Professional L Headband

MOUNTING OPTION

In a lightweight professional design for fatigue-free work

2



HR 2.5x/340mm



HR 2.5x/420mm



HR 2.5x/520mm



HRP 3.5x/420mm



HRP 4x/340mm



HRP 6x/340mm

MAGNIFICATION AND WORKING DISTANCE

Choose the smallest magnification that fulfils your requirements. The lower the magnification, the easier it is to use the loupes due to the larger field of view and better depth of focus. The working distance depends on personal preference.

3



ILLUMINATION

Choose between the LoupeLight2 and the ML4 LED.



LoupeLight2
for S-FRAME

Compact and lightweight.
Coaxial illumination.



ML4 LED HeadLight
for Professional L Headband

Coaxial illumination.
Flexible power source options.

4

POWER SOURCE



mPack mini
Portable
power pack



mPack
Portable
power pack



mPack UNPLUGGED
Wireless.
Headband-mounted